



Opportunities, challenges, and economic drivers for start-to-finish continuous biomanufacturing

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY15169-13Jul20-PT



Opportunities, challenges, and economic drivers for start-to-finish continuous biomanufacturing

Avril Vermunt, Andreas Castan, Oliver Hardick, Magnus Bergman, Hans Blom, Ellen Johnson, Lotta Molander, and Günter Jagschies. GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden

Introduction

Biomufacturing today is informed by our past and considers our current state

Continuous unit operations have been used successfully in manufacturing of several approved products with a collective revenue of \$25 billion annually. Many large biopharmaceutical organizations are exploring potential solutions to drive continuous operations on the manufacturing floor. Continuous biomanufacturing carries the hope for achieving the next level of efficiency through significant increases in productivity, floor space reduction, and lower capital cost for manufacturing facilities as well as lower overall production costs. Combining many innovative, promising unit operation technologies with single-use equipment can offer significant economic advantages over batch processing in traditional stainless steel installations.

Adopting advanced process design will need an accurate business case and risk profile for implementation. Within the risk profile, assurance of patient safety and supply through product quality demonstration will be critical. Optimizing the operating space of an entire process, versus a unit operation, will call for thoughtful automation and control strategies. Simplification of equipment design and operator interface points will also be key to success.

This poster will highlight trade-offs and benefits associated with the adoptions of connected and continuous manufacturing.

One of the first challenges to recognize is that decision making and cost analyses combine multiple sources of savings (Fig 1), and often complicate the view of each one or combination of them.




Fig 1. Example sources of savings.

Benefits and challenges

Examples will be highlighted for the following continuous and connected biomanufacturing steps:

- High cell density WCBs
- N-1 and N perfusion
- Protein A fibro units
- Continuous chromatography
- Automation and controls
- CAPEX strategies

Cell culture

Modernize processes by implementing high cell density WCBs and intensify N-1 through perfusion (1)




Fig 2. Reduction of time visualized high cell density WCB and perfusion N-1.

There are benefits related to operations time reduction (Fig 2), but there are challenges to processing the increased output. How can continuous biomanufacturing address this?

Continuous unit operations help offset cost implication of adopting modern manufacturing practices

Assumptions:

- 6 x 2000 L at 10 g/L and 3 x 500 L at 6 g/L/d, adjusted for similar annual product output of ~ 1340 kg mAb product
- 1.1 reactor volumes of media in batch mode, two reactor volumes of media per day in perfusion, 10% bleed
- Similar culture media composition and cost/L

The cost of increased titer or output can be balanced with the savings of a continuous operating mode (Fig 3) as long as media costs are managed (2) (Fig 4).




Fig 3. Number of single-use bioreactor bags for batch process vs perfusion process. Fig 4. Cost of media per gram mAb for batch vs perfusion process.

mAb capture

Protein A fibro units can operate in batch mode connected to upstream and purify the product within 24 h

Immediate mass transfer enables high productivity (Fig 5 and 6). By operating Protein A fibro units in batch mode connected to upstream you can purify the product within 24 h.




Fig 5. Dynamic binding capacity and pressure drop. Fig 6. An industrial collaborator achieved 460 g/L/h productivity purifying 40 g mAb with equivalent product quality with a 10 mL unit in one day.

Process productivity and flexibility improved because:

- Process performance is independent of flow rate and tolerant of process variation.
- Technology supports connected manufacturing and scale-out.

Continuous chromatography

For some scales and campaign approaches, continuous chromatography offers value in purification resin savings (Table 1).

Table 1. Manufacturing scenarios with process modes and example resin costs

Manufacturing scenarios	One-time campaigns 500 g to 10 kg	Medium length campaigns 10 to 500 kg/yr	All-day/all-year operations 200 to 10 000 kg/yr
Batch mode	Resin cost: depends on scale and output	Resin cost: \$6-8/g Define process for best capacity and maximum resin utilization.	Resin cost: \$3-6/g
Batch + continuous modes	-	Resin cost: ~\$3/g Optimize number of runs with respect to best capacity resin	Resin cost: \$2-3/g
Continuous mode	-	Resin cost: \$2-3/g Overcome time-based lifetime (or run/cycle-based lifetime) constraint with PCC	Resin cost: ~\$2/g Optimize capacity with PCC

Acronyms

API active pharmaceutical ingredient	OPEX operating expenditure
CAPEX capital expenditure	PAT process analytical technology
EMA European Medicines Agency	PCC periodic counter-current chromatography
FDA the US Food and Drug Administration	PD process development
GMP good manufacturing practice	ROI return on investment
HVAC heating, ventilation, and air conditioning	WCD working cell bank

Operations

Flexible integration

Flexible integration approaches based on automation expertise and bioprocess knowledge accomplish GMP control for continuous and connection biomanufacturing (Fig 7 and Table 2).

Table 2. Levels of automation integration

- Unit level**
 - Islands of automation for single unit operations
 - Real-time monitoring, tracking, and optimization
 - Reduced variability
 - GMP documentation
- Process level**
 - Integration at process level for central control
 - Manufacturing consistency
 - Turn data insights into outcomes
- Plant level**
 - Distributed control system (DCS) for facility wide control
 - Site data management
 - Process, utilities, and HVAC



Fig 7. Visual examples of integration levels.

Reduced CAPEX

Reduced CAPEX is a culmination of modernization, intensification, and industrialization efforts, and is enabled by continuous and connected biomanufacturing (Fig 8).




Fig 8. CAPEX as published by biopharma companies and CMOs (blue), recently constructed/announced single-use or hybrid facilities (orange). 15 plus years of facility construction history, corrected for inflation, typically showing greenfield construction

References

- Application note: One step seed train from one cell bank vial to 2000 L production bioreactor, GE Healthcare, 29160932, Edition AB
- Xu, S. Bioreactor productivity and media cost comparison for different intensified cell culture processes, *Biotechnol. Proc.* **33**(4), 867-878 (2017)
- With permission: Jagschies, G. et al., *Biopharmaceutical Processing: Development, Design, and Implementation of Manufacturing Processes*, Elsevier (2018) (graph updated)

Summary

Benefits with start-to-finish continuous biomanufacturing when done right

Indirect financial benefits:

- Revenue and pipeline growth due to flexibility of scale and product mix
- Commercial scale for manufacturing drives the clinical scale
- Reduction of analytical complexity and cost
- Reduction of labor cost

Direct financial benefits:

- Smaller and less complex equipment
- Increased facility productivity
- Lower direct cost for supplied materials and consumables
- Lower cost of use for supplied equipment, consumables, and materials

Implementing start-to-finish processes requires a balanced approach

Table 3. Pros and cons for start-to-finish processing

Pros	Cons
Right size scale and CAPEX reduction	OPEX minimally lowered
Production scale drives CAPEX reduction	SU downstream adds risk with minimal improvement to ROI
Necessitates better operating space definition	Cell culture media cost requires management
Enables varied product mix	Overall risk is increased
Reduction of non-value added activities	Resource shift for development, tech transfer, and validation

Biomufacturing of the future

Adoption of connected and continuous processing has been slow as basic changes in manufacturing paradigms take several years, in part due to inherent industry conservatism and the high regulation of the industry. The top concerns include contamination risk, process control challenge, and operational complexity. Continuous and connected processing is viewed as complex and susceptible to problems. These perceptions are not in line with the advances in the industry.

Both FDA and EMA have indicated openness due to the view of better process control, product quality, and easier application of PAT. Upcoming start to finish processes can understanding the challenges, as well as benefits, and balancing decisions to ensure successful results.